

## samaritan® PAD Family







	SAM 350P with CPR Coaching	SAM 450P with CPR Rate Advisor™
Integrated CPR Rate Advisor™ ICG-based feedback on CPR rate		<b>✓</b>
CPR coaching with metronome	<b>✓</b>	<b>✓</b>
Audible/visual prompts to maximize CPR hands-on time		<b>✓</b>
Real-time audible CPR rate feedback		<b>✓</b>
Real-time visual CPR rate feedback		<b>✓</b>
Audible and visual indicators for ease of use in a noisy environment	<b>✓</b>	<b>✓</b>
Compact, portable, and lightweight 1.1 kg/2.4 lbs	<b>✓</b>	<b>✓</b>
Highest level of protection from dust and water ingress in the industry (IP56)	<b>✓</b>	<b>✓</b>
Integrated battery and electrodes for one-change maintenance	<b>✓</b>	<b>✓</b>
Four-year electrode and battery life	<b>✓</b>	<b>✓</b>
Proprietary low escalating energy protocol (SCOPE™)	<b>✓</b>	<b>✓</b>
Semi-automatic shock delivery	<b>✓</b>	<b>✓</b>
Pediatric compatible	<b>✓</b>	<b>✓</b> *
90-minute ECG recording time	<b>✓</b>	<b>✓</b>
Automatic self test	<b>✓</b>	<b>✓</b>
Field upgradeable software to latest AHA/ERC guidelines	<b>✓</b>	V
Carry case included		<b>✓</b>

<sup>\*</sup> CPR Rate Advisor is deactivated when the Pediatric-Pak is in use.

## EMEA/APAC

HeartSine Technologies, Ltd. 203 Airport Road West Belfast, Northern Ireland BT3 9ED Tel: +44 28 9093 9400 Fax: +44 28 9093 9401

info@heartsine.com



## **U.S./Americas**

HeartSine Technologies, Inc. 121 Friends Lane, Suite 400 Newtown, PA 18940 Toll Free: (866) 478 7463 Tel: +1 215 860 8100 Fax: +1 215 860 8192 info@heartsine.com

The samaritan PAD 350P described in this brochure meets the applicable European Medical Directive requirements.



CUL)us UL Classified. See complete marking on product.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a licensed practitioner.

© 2015 HeartSine Technologies, Inc. All rights reserved. Hoo9-032-336-0

**Lifesaving,** Pure and Simple

www.heartsine.com







